

## Product datasheet for **RC217574L1V**

### ADAMTS15 (NM\_139055) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ADAMTS15 (NM_139055) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ADAMTS15
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_139055
ORF Size:	2850 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217574).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_139055.1</a>
RefSeq Size:	2853 bp
RefSeq ORF:	2853 bp
Locus ID:	170689
UniProt ID:	<a href="#">Q8TE58</a>
Cytogenetics:	11q24.3
Protein Families:	Druggable Genome, Protease, Secreted Protein
MW:	101.5 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the ADAMTS (a disintegrin and metalloproteinase with thrombospondin motifs) protein family. ADAMTS family members share several distinct protein modules, including a propeptide region, a metalloproteinase domain, a disintegrin-like domain, and a thrombospondin type 1 (TS) motif. Individual members of this family differ in the number of C-terminal TS motifs, and some have unique C-terminal domains. The encoded preproprotein is proteolytically processed to generate the mature enzyme, which may play a role in versican processing during skeletal muscle development. This gene may function as a tumor suppressor in colorectal and breast cancers. [provided by RefSeq, May 2016]