

Product datasheet for **RC217575L3V**

LOXL3 (NM_032603) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	LOXL3 (NM_032603) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LOXL3
Synonyms:	LOXL
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032603
ORF Size:	2259 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217575).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032603.2
RefSeq Size:	3121 bp
RefSeq ORF:	2262 bp
Locus ID:	84695
UniProt ID:	P58215
Cytogenetics:	2p13.1
Protein Families:	Druggable Genome
MW:	80.3 kDa


[View online »](#)

Gene Summary:

This gene encodes a lysyl oxidase, which likely functions as an amine oxidase and plays a role in the formation of crosslinks in collagens and elastin. Deletion of the related gene in mouse causes neonatal mortality with cleft palate, spine deformity, and defects in collagen organization. A mutation in this gene was found in a family with Stickler syndrome. [provided by RefSeq, Sep 2016]