

Product datasheet for **RC217579L3V**

CDR1 (NM_004065) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CDR1 (NM_004065) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CDR1
Synonyms:	CDR; CDR34; CDR62A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_004065
ORF Size:	786 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217579).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004065.2 , NP_004056.2
RefSeq Size:	1299 bp
RefSeq ORF:	789 bp
Locus ID:	1038
UniProt ID:	P51861
Cytogenetics:	Xq27.1
MW:	31.1 kDa



[View online »](#)

Gene Summary:

Autoantibodies directed against the protein encoded by this intronless gene have been found in some patients with paraneoplastic cerebellar degeneration. The encoded protein contains several hexapeptide repeats. [provided by RefSeq, Jan 2010]