

Product datasheet for RC217580L4V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

FOXP3 (NM_014009) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: FOXP3 (NM_014009) Human Tagged ORF Clone Lentiviral Particle

Symbol: FOXP3

Synonyms: AIID; DIETER; IPEX; JM2; PIDX; XPID

Mammalian Cell

Selection:

Puromycin

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

Tag: mGFP

ACCN: NM_014009 **ORF Size:** 1293 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC217580).

OTI Disclaimer:

Sequence:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeg: NM 014009.2

 RefSeq Size:
 1869 bp

 RefSeq ORF:
 1296 bp

 Locus ID:
 50943

 UniProt ID:
 Q9BZS1

 Cytogenetics:
 Xp11.23

Protein Families: Transcription Factors

MW: 47.1 kDa







Gene Summary:

The protein encoded by this gene is a member of the forkhead/winged-helix family of transcriptional regulators. Defects in this gene are the cause of immunodeficiency polyendocrinopathy, enteropathy, X-linked syndrome (IPEX), also known as X-linked autoimmunity-immunodeficiency syndrome. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]

Product images:



[RC217580L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC217580L4V particle to overexpress human FOXP3-mGFP fusion protein.



[RC217580L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC217580L4V particle to overexpress human FOXP3-mGFP fusion protein.