

Product datasheet for RC217587L4

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Argininosuccinate Lyase (ASL) (NM_001024944) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Argininosuccinate Lyase (ASL) (NM_001024944) Human Tagged Lenti ORF Clone

Tag: mGFP

Symbol: Argininosuccinate Lyase

Synonyms: ASAL

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this of


Sequence:

The ORF insert of this clone is exactly the same as(RC217587).

Restriction Sites: Sgfl-Rsrll

Cloning Scheme:





 $[\]ensuremath{^*}$ The last codon before the Stop codon of the ORF.

ACCN: NM_001024944

ORF Size: 1332 bp



Argininosuccinate Lyase (ASL) (NM_001024944) Human Tagged Lenti ORF Clone - RC217587L4

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 001024944.1</u>

RefSeq Size: 2001 bp RefSeq ORF: 1335 bp Locus ID: 435

 UniProt ID:
 P04424

 Cytogenetics:
 7q11.21

Protein Pathways: Alanine, aspartate and glutamate metabolism, Arginine and proline metabolism, Metabolic

pathways

MW: 49.3 kDa


Gene Summary: This gene encodes a member of the lyase 1 family. The encoded protein forms a cytosolic

homotetramer and primarily catalyzes the reversible hydrolytic cleavage of argininosuccinate into arginine and fumarate, an essential step in the liver in detoxifying ammonia via the urea cycle. Mutations in this gene result in the autosomal recessive disorder argininosuccinic aciduria, or argininosuccinic acid lyase deficiency. A nontranscribed pseudogene is also located on the long arm of chromosome 22. Alternatively spliced transcript variants encoding

different isoforms have been described. [provided by RefSeq, Jul 2008]



Product images:



Circular map for RC217587L4