

Product datasheet for **RC217610L3V**

KIF21A (NM_017641) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	KIF21A (NM_017641) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KIF21A
Synonyms:	CFEOM1; FEOM1; FEOM3A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017641
ORF Size:	4983 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217610).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017641.2
RefSeq Size:	6335 bp
RefSeq ORF:	4986 bp
Locus ID:	55605
UniProt ID:	Q7Z4S6
Cytogenetics:	12q12
Protein Families:	Druggable Genome
MW:	185.3 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the KIF4 subfamily of kinesin-like motor proteins. The encoded protein is characterized by an N-terminal motor domain a coiled-coil stalk domain and a C-terminal WD-40 repeat domain. This protein may be involved in microtubule dependent transport. Mutations in this gene are the cause of congenital fibrosis of extraocular muscles-1. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Mar 2010]