

Product datasheet for **RC217616L4V**

BRD3 (NM_007371) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	BRD3 (NM_007371) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BRD3
Synonyms:	ORFX; RING3L
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007371
ORF Size:	2178 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217616).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007371.3
RefSeq Size:	5673 bp
RefSeq ORF:	2181 bp
Locus ID:	8019
UniProt ID:	Q15059
Cytogenetics:	9q34.2
Domains:	BROMO
Protein Families:	Protein Kinase



[View online »](#)

MW: 80 kDa

Gene Summary: This gene was identified based on its homology to the gene encoding the RING3 protein, a serine/threonine kinase. The gene localizes to 9q34, a region which contains several major histocompatibility complex (MHC) genes. The function of the encoded protein is not known. [provided by RefSeq, Jul 2008]