

Product datasheet for **RC217617L3V**

DOK1 (NM_001381) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DOK1 (NM_001381) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DOK1
Synonyms:	P62DOK; pp62
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001381
ORF Size:	1443 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217617).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001381.2
RefSeq Size:	1908 bp
RefSeq ORF:	1446 bp
Locus ID:	1796
UniProt ID:	Q99704
Cytogenetics:	2p13.1
Domains:	PH, IRS
Protein Families:	Druggable Genome



[View online »](#)

MW: 52.2 kDa

Gene Summary: The protein encoded by this gene is part of a signal transduction pathway downstream of receptor tyrosine kinases. The encoded protein is a scaffold protein that helps form a platform for the assembly of multiprotein signaling complexes. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2016]