

Product datasheet for **RC217620L4V**

SLC5A8 (NM_145913) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	SLC5A8
Synonyms:	AIT; SMCT; SMCTI
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145913
ORF Size:	1830 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC217620).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_145913.2
RefSeq Size:	3286 bp
RefSeq ORF:	1833 bp
Locus ID:	160728
UniProt ID:	Q8N695
Cytogenetics:	12q23.1-q23.2



Domains: SSF

Protein Families: Transmembrane

MW: 66.6 kDa

Gene Summary: SLC5A8 has been shown to transport iodide by a passive mechanism (Rodriguez et al., 2002 [PubMed 12107270]) and monocarboxylates and short-chain fatty acids by a sodium-coupled mechanism (Gopal et al., 2004 [PubMed 15322102]). In kidney, SLC5A8 functions as a high-affinity sodium-coupled lactate transporter involved in reabsorption of lactate and maintenance of blood lactate levels (Thangaraju et al., 2006 [PubMed 16873376]).[supplied by OMIM, Dec 2008]