

Product datasheet for **RC217625**

BCNP (FAM129C) (NM_001098524) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BCNP (FAM129C) (NM_001098524) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	BCNP
Synonyms:	BCNP1; FAM129C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC217625 representing NM_001098524
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGGCCTGACCGGAAGGAGGTGCCCTGAGCCGAGGAACGAGCGGTGGTCTGGGGAAGGAAGAG
GAGCCCCGGGAGACGACAGCAGCATGGGTGGGCGCCTTCGAGCCCTCTGGACAAGCAGCAGCGGCAGCA
CCTAAGGGGTCAGGTGGACACCCTGCTGAGGAACTTCTGCCTTGTACCGTGGGCACTGGCAGCGTCT
GTCCTGCGGCAGATCTCTCGAGAGCTGGGCCCTCAGGAGCCGACCGGAAGCCAGTTGCTACGCAGCAAAA
AGCTGCCCGAGTCCGTGAGCACCGAGGACCCCTGACCCAGCTTCGGGGCCACCCACCCCGGTGGCAGCC
GATCTTCTGTGTTCTGCGTGGGACGGCCCTAGAGTGGTTCAGCCACAAGGAGGAATATGAAAACGGG
GGCCACTGCCTGGCTCAACAGCCCTGACAGGATACACGCTCCTGACTTCCAGCGAGAATATCTCCGCC
TTTTGGATGCTCTGCCCTGAATCCTTGGGAGACCATACTCAGGAAGAGCCTGACTCCCTTTGGAAGT
GCCTGTGAGCTTCCCGCTGTTCTGCAGCACCCCTTCCGCCGGCACCTCTGCTTCTGCAGCCACCAGG
GAGGCACAGCATGCCTGGAGGCTGGCCCTGCAGGGTGGCATCCGGCTTCAGGGCAGATCCTGCAGCGAA
GCCAGGCCCTGCTGCCCGGCCCTTCTGGACGCCGTCGACTCTACCGGCAGCACCAAGGCCACTTTGG
CGACGACGACGTGACCCTAGGCTCAGACGCCGAGGTGCTGACCGCGGTGCTGATGCGGGAGCAACTTCCC
GCGCTGCGAGCCAGACCCTTCTGGCCTGCAGGGGGCAGGCCGCGCCCGCCTGGGCCTGGACCAGC
TTCTAGACGCCGTTACGCAGCTGTCTGGCCGGGGCCTCCGCCGGGCTCTGCGCCTTCCAGCCGAAAA
GGACGAGCTGCTGCGTGCCTGGAGAAGACGATCCGCCGGACGTGGACCAGTCTGCGGCAGCGGGCG
CGTGTGGCGGGGGCTGAGGACGGATACAGGGGACCGCTCAGTCTGCTGCGCCGGGAGGTGGACC
CGCAGCTGCCCGGGTCTGCAGACCCTGCTGCCACCGTGGAGCCCTGCTCGAGCGGTGCGGACCCCT
CCTGGCTCAAGGCATGGACCGACTGTCCACCGCTGCGCCAGAGCCCTCAGGCACGCGGCTGCGCAGG
GAGGTTTACTCATTTGGGAGATGCCGTGGGACTTGGCGCTGATGCAGACATGCTACCGTGAAGCCGAGC
GGAGCCGGGGCGCTTGGGGCAGCTGGCAGCACCGTTTGGCTTTCTGGGATGCAGAGCCTCGTGTTTGG
GGCCCAAGATCTTGACAGCAGCTCATGGCTGACGCCGTGGCCACCTTCTGCAGCTGGCTGACCAGTGT
CTGACGACGGCCCTCAACTGTGACCAGGCTGCCAGAGGCTGGAGAGAGTCAGGGGGCGCGTGTGAAGA
AATTCAAATCGGACAGCGGGTTGGCGCAGAGGAGGTTCCAGGCTGGGGTCTCTGCATCTTTTTACC
TTTTGTGCTGAGCCAACCGAGCCAGGCTGCAAAAAGGAGCTGCCTGAGTTCAGGGGGATGCTTGGC
GTGGGCAGCCAGGCTCTGACCACTGAGGGCATCTATGAGGACGTATCCGGGGTCTTGTGTCAGAGGA
TTGACCAAGAATTGAAAAGACCCTTGGTGCCAATGATGTATCCTGCACTCTGGACGGCTGCTTGGAGGT
CCCATGGGAACAGGAGGGAGCAGATGAGGAACTGAGGCTGAGCGGGAAGGAGGGGCTTGTCCCAGGCAG
CCAGACTCTGGTCCAGATCCAGCCACTCTGCCACCGCCTTCTCCAGGAACATTCCGGAGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC217625 representing NM_001098524
Red=Cloning site Green=Tags(s)


MGPDRKEVPLSRGTQAVVVGKGRGAPGDDSSMGGRPSSPLDKQQRQHLRGQVDTLLRNFLPCYRGQLAAS
 VLRQISRELGPQEPTGSQLLRSKLLPRVREHRGPLTQLRGHPRWQPIFCVLRGDGRLEWFSHKEEYENG
 GHCLGSTALTGYLLTSQREYLRLLDALCPESLGDHTQEEPDSLLEVPVSPFLFLQHPFRRHLCSAATR
 EAQHAWRLALQGGIRLQGTVLQRSQAPAAAFLLDAVRLYRQHGHFGDDDDVTLGSDAEVL TAVLMREQLP
 ALRAQTLPLGLRGAGRARAWATELLDAVHA AVLAGASAGLCAFQPEKDELLASLEKTIRPDVDQLLRQRA
 RVAGRLRTDIRGPLESCLRREVDPQLPRVVQTLRLTVEASLEAVRLLAQGMDRLSHRLRQSPSGTRLRR
 EYVSFGEMPWDLALMQTCYREAERSRGLGQLAAPFGFLGMQSLVFQAQDLAQQLMADAVATFLQLADQC
 LTTALNCDQAAQRLERVRGRVLKFKSDSGLAQRRFIRGWGLCIFLPFVLSQLEPGCKKELPEFEGDVLA
 VGSQALTEGIYEDVIRGCLLQRIDQELKKT LGANDV SCTL DGCLEVPWEQEGADEETEAREGGACPRQ
 PDSGAQIQPLCPPSPGTFRS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8071_e02.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001098524

ORF Size: 1953 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001098524.1](#), [NP_001091994.1](#)

RefSeq Size: 2233 bp

RefSeq ORF: 1956 bp


Locus ID: 199786

UniProt ID: [Q86XR2](#)

Cytogenetics: 19p13.11

MW: 72 kDa

Product images:



Circular map for RC217625