

## Product datasheet for **RC217648A1V**

### Human Sclerostin (SOST) (NM\_025237) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human Sclerostin (SOST) (NM\_025237) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Sclerostin  
**Synonyms:** CDD; DAND6; SOST1; VBCH  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC217648 representing NM\_025237  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGCC**

ATGCAGCTCCCACTGGCCCTGTGTCTCGTCTGCCTGCTGGTACACACAGCCTTCCGTGTAGTGGAGGGCC  
AGGGGTGGCAGGCGTTCAAGAATGATGCCACGGAATCATCCCGAGCTCGGAGAGTACCCCGAGCCTCC  
ACCGGAGCTGGAGAACAAAGACCATGAACCGGGCGGAGAACGGAGGGCGGCCCTCCCAACACCCCTTT  
GAGACCAAAGACGTGTCCGAGTACAGCTGCCCGGAGCTGCACTTACCCGCTACGTGACCGATGGGCCGT  
GCCGACGCGCAAGCCGGTCACCGAGCTGGTGTCTCCGGCCAGTGCAGGCCCGGCCGCGCTGCTGCCAA  
GCCCATCGGCCCGGCAAGTGGTGGCGACCTAGTGGGCCGACTTCCGCTGCATCCCGACCGCTACCGC  
GCGCAGCGCGTGCAGCTGCTGTGTCCCGGTGGTGGGCGCCGCGCGCGCAAGGTGCGCCTGGTGGCCT  
CGTGCAAGTGCAAGCGCCTACCCGCTTCCACAACCAAGTCGGAGCTCAAGGACTTCGGGACCGAGGCCGC  
TCGGCCGCGAGAAGGGCCGAAGCCGCGGCCCGCGCCCGGAGCGCAAGCCAACCGAGCCGAGCTGGAG  
AACGCCTAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAATGATATCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC217648 representing NM\_025237  
Red=Cloning site Green=Tags(s)

MLPLALCLVCLLVHTAFRVVEGQGWQAFKNDATETIPELGEYPEPPPELENNKTMNRAENGRRPPHPF  
ETKDVSEYSCRELHFTRYVTDGPCRSAPVTELVCSGQCGPARLLPNAIGRGKWWRPSGPDFRCIPDRYR  
AQRVQLLCPGGEAPRARKVRLVASCKCKRLTRFHNQSELKDFGTEAARPQKGRKPRPRARSAKANQAELE  
NAY

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

|                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Species:</b>      | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Serotype:</b>     | AAV-2                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>ACCN:</b>         | NM_025237                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>ORF Size:</b>     | 639 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Buffer:</b>       | PBS with 0.001% Pluronic F68                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Stability:</b>    | AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus. |
| <b>RefSeq:</b>       | <a href="#">NM_025237.2</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq Size:</b>  | 2322 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>RefSeq ORF:</b>   | 642 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Locus ID:</b>     | 50964                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>UniProt ID:</b>   | <a href="#">Q9BQB4</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Cytogenetics:</b> | 17q21.31                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>MW:</b>           | 21.5 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |