

Product datasheet for RC217666

OR52M1 (NM_001004137) Human Tagged ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | OR52M1 (NM_001004137) Human Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | OR52M1 |
| Synonyms: | OR11-11; OR52M1P; OR52M3P |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| ORF Nucleotide Sequence: | >RC217666 representing NM_001004137 Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTCACTTTTCATAATGTCTGCTCAGTACCCAGCTCCTTCTGGCTCACTGGCATCCCAGGGCTGGAGT
CCCTACACGTCTGGCTCTCCATCCCCTTTGGCTCCATGTACCTGGTGGCTGTGGTGGGAATGTGACCAT
CCTGGCTGTGGTAAAGATAGAACGCAGCCTGCACCAGCCCATGTACTTTTCTTGTGCATGTTGGCTGCC
ATTGACCTGGTTCTGTCTACTTCCACTATACCCAACTTCTGGGAATCTTCTGGTTCGGTGCTTGTGACA
TTGGCCTGGAGCCTGCTTGGGCCAAATGTTCTTATCCACTGCTTGGCCACTGTTGAGTCAGGCATCTT
CCTTGCCATGGCTTTTGATCGCTACGTGGCCATCTGCAACCCACTACGTCATAGCATGGTGTCACTTAT
ACAGTGGTGGGTCGTTTGGGGCTTGTCTCTCCTCCGGGTGTTCTCTACATTGGACCTCTGCCTCTGA
TGATCCGCTGCGGCTGCCCTTTATAAAACCCATGTTATCTCCCACTCCTACTGTGAGCACATGGCTGT
AGTTGCCCTTGACATGTGGCGACAGCAGGGTCAATAATGTCTATGGGCTGAGCATCGGCTTTCTGGTGTG
ATCCTGGACTCAGTGGCTATTGCTGCATCCTATGTGATGATTTTCAGGGCCGTGATGGGGTTAGCCACTC
CTGAGGCTAGGCTTAAACCCCTGGGGACATGCGCTTCTCACCTCTGTGCCATCCTGATCTTTTATGTTCC
CATTGCTGTTTCTCCCTGATTCACCGATTTGGTCAGTGTGTGCCTCCTCCAGTCCACACTCTGCTGGCC
AACTTCTATCTCCTCATTCTCCAATCCTCAATCCCATTGTCTATGCTGTTTCGACCAAGCAGATCCGAG
AGAGCCTTCTCCAAATACCAAGGATAGAATGAAGATTAGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC217666 representing NM_001004137
 Red=Cloning site Green=Tags(s)

MLTFHNVCSVPSSFWLTGIPGLES�HVLSIPFGSMYLVAVVGNVTILAVVKIERSLHQPMYFFLCMLAA
 IDLVLSTSTIPKLLGIFWFGACDIDLDAQLGQMFLIHCFATVESGIFLAMAFDRYVAICNPLRHSMVLT
 TVVGRLLVSLLRGVLYIGPLPLMIRLRLPLKTHVISHSYCEHMAVVALTCGDSRVNNVYGLSIGFLVL
 ILDSVAIAASYVMIFRAVMGLATPEARLKTGTGCASHLCAILIFYVPIAVSSLIHRFGQCVPPPVTLLA
 NFYLLIPPILNPIVYAVRTKQIRESLLQIPRIEMKIR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8001_h03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001004137

ORF Size: 951 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001004137.1](#), [NP_001004137.1](#)

RefSeq Size: 954 bp

RefSeq ORF: 954 bp

Locus ID: 119772

UniProt ID: [Q8NGK5](#)

Cytogenetics: 11p15.4


Protein Families: Transmembrane

Protein Pathways: Olfactory transduction

MW: 34.9 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC217666