

Product datasheet for **RC217712L3V**

HMGCLL1 (NM_019036) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	HMGCLL1
Synonyms:	bA4I8P12.1; er-cHL
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_019036
ORF Size:	1110 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC217712).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_019036.2 , NP_061909.2
RefSeq Size:	2571 bp
RefSeq ORF:	1113 bp
Locus ID:	54511
UniProt ID:	Q8TB92
Cytogenetics:	6p12.1



MW: 39.3 kDa

Gene Summary: Non-mitochondrial 3-hydroxymethyl-3-methylglutaryl-CoA lyase that catalyzes a cation-dependent cleavage of (S)-3-hydroxy-3-methylglutaryl-CoA into acetyl-CoA and acetoacetate, a key step in ketogenesis, the products of which support energy production in nonhepatic animal tissues.[UniProtKB/Swiss-Prot Function]