

## Product datasheet for **RC217729L1V**

### PHC3 (NM\_024947) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PHC3 (NM_024947) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PHC3
Synonyms:	EDR3; HPH3
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_024947
ORF Size:	2949 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217729).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024947.3</a>
RefSeq Size:	12687 bp
RefSeq ORF:	2988 bp
Locus ID:	80012
UniProt ID:	<a href="#">Q8NDX5</a>
Cytogenetics:	3q26.2
Domains:	SAM
MW:	106.2 kDa



[View online »](#)

**Gene Summary:**

Component of a Polycomb group (PcG) multiprotein PRC1-like complex, a complex class required to maintain the transcriptionally repressive state of many genes, including Hox genes, throughout development. PcG PRC1 complex acts via chromatin remodeling and modification of histones; it mediates monoubiquitination of histone H2A 'Lys-119', rendering chromatin heritably changed in its expressibility.[UniProtKB/Swiss-Prot Function]