

## Product datasheet for **RC217730L4V**

### KIF25 (NM\_005355) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	KIF25 (NM_005355) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KIF25
Synonyms:	KNLSL3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_005355
ORF Size:	996 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217730).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_005355.2</a> , <a href="#">NP_005346.3</a>
RefSeq Size:	1354 bp
RefSeq ORF:	999 bp
Locus ID:	3834
UniProt ID:	<a href="#">Q9UIL4</a>
Cytogenetics:	6q27
Protein Families:	Druggable Genome
MW:	35.1 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a member of the kinesin-like protein family. Protein family members are microtubule-dependent molecular motors that transport organelles within cells and move chromosomes during cell division. However, the particular function of this gene product has not yet been determined. Two alternatively spliced transcript variants which encode products have been described. Other splice variants have been found that lack exon 2 and the initiation codon for translation. [provided by RefSeq, Jul 2008]