

Product datasheet for **RC217739L3V**

PER3 (NM_016831) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PER3 (NM_016831) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PER3
Synonyms:	FASPS3; GIG13
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_016831
ORF Size:	3603 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217739).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016831.1
RefSeq Size:	6203 bp
RefSeq ORF:	3606 bp
Locus ID:	8863
UniProt ID:	P56645
Cytogenetics:	1p36.23
Domains:	PAS
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

Protein Pathways: Circadian rhythm - mammal

MW: 131.7 kDa

Gene Summary: This gene is a member of the Period family of genes and is expressed in a circadian pattern in the suprachiasmatic nucleus, the primary circadian pacemaker in the mammalian brain. Genes in this family encode components of the circadian rhythms of locomotor activity, metabolism, and behavior. This gene is upregulated by CLOCK/ARNTL heterodimers but then represses this upregulation in a feedback loop using PER/CRY heterodimers to interact with CLOCK/ARNTL. Polymorphisms in this gene have been linked to sleep disorders. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2014]