

Product datasheet for **RC217744L3V**

CARHSP1 (NM_001042476) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CARHSP1 (NM_001042476) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CARHSP1
Synonyms:	CRHSP-24; CRHSP24; CSDC1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001042476
ORF Size:	441 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217744).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001042476.1 , NP_001035941.1
RefSeq Size:	3029 bp
RefSeq ORF:	444 bp
Locus ID:	23589
UniProt ID:	Q9Y2V2
Cytogenetics:	16p13.2
Protein Families:	Transcription Factors
MW:	15.7 kDa



[View online »](#)

Gene Summary:

Binds mRNA and regulates the stability of target mRNA. Binds single-stranded DNA (in vitro).
[UniProtKB/Swiss-Prot Function]