

## Product datasheet for **RC217756L4V**

### BRAP (NM\_006768) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	BRAP (NM_006768) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BRAP
Synonyms:	BRAP2; IMP; RNF52
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006768
ORF Size:	1776 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217756).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006768.2</a>
RefSeq Size:	2063 bp
RefSeq ORF:	1779 bp
Locus ID:	8315
UniProt ID:	<a href="#">Q7Z569</a>
Cytogenetics:	12q24.12
Protein Families:	Druggable Genome
MW:	67.1 kDa


[View online »](#)

**Gene Summary:**

The protein encoded by this gene was identified by its ability to bind to the nuclear localization signal of BRCA1 and other proteins. It is a cytoplasmic protein which may regulate nuclear targeting by retaining proteins with a nuclear localization signal in the cytoplasm. [provided by RefSeq, Jul 2008]