

## Product datasheet for **RC217780L4V**

### TDRKH (NM\_001083964) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TDRKH (NM_001083964) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TDRKH
Synonyms:	TDRD2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001083964
ORF Size:	1548 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217780).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001083964.1</a> , <a href="#">NP_001077433.1</a>
RefSeq Size:	2719 bp
RefSeq ORF:	1551 bp
Locus ID:	11022
UniProt ID:	<a href="#">Q9Y2W6</a>
Cytogenetics:	1q21.3
Protein Families:	Transmembrane
MW:	57 kDa



[View online »](#)

**Gene Summary:**

Participates in the primary piRNA biogenesis pathway and is required during spermatogenesis to repress transposable elements and prevent their mobilization, which is essential for the germline integrity. The piRNA metabolic process mediates the repression of transposable elements during meiosis by forming complexes composed of piRNAs and Piwi proteins and govern the methylation and subsequent repression of transposons. Required for the final steps of primary piRNA biogenesis by participating in the processing of 31-37 nt intermediates into mature piRNAs. May act in pi-bodies and piP-bodies by transferring piRNA precursors or intermediates to or between these granules.[UniProtKB/Swiss-Prot Function]