

Product datasheet for **RC217810L4V**

ENDOD1 (NM_015036) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ENDOD1 (NM_015036) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ENDOD1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_015036
ORF Size:	1500 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217810).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015036.1
RefSeq Size:	4650 bp
RefSeq ORF:	1503 bp
Locus ID:	23052
UniProt ID:	O94919
Cytogenetics:	11q21
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Apoptosis
MW:	54.8 kDa



[View online »](#)

Gene Summary: May act as a DNase and a RNase.[UniProtKB/Swiss-Prot Function]