

## Product datasheet for **RC217815L3V**

### STARD13 (NM\_052851) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	STARD13 (NM_052851) Human Tagged ORF Clone Lentiviral Particle
Symbol:	STARD13
Synonyms:	ARHGAP37; DLC2; GT650; LINC00464
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_052851
ORF Size:	2985 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217815).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_052851.2</a>
RefSeq Size:	5866 bp
RefSeq ORF:	2988 bp
Locus ID:	90627
UniProt ID:	<a href="#">Q9Y3M8</a>
Cytogenetics:	13q13.1-q13.2
Domains:	RhoGAP, START
MW:	111.2 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a protein which contains an N-terminal sterile alpha motif (SAM) for protein-protein interactions, followed by an ATP/GTP-binding motif, a GTPase-activating protein (GAP) domain, and a C-terminal STAR-related lipid transfer (START) domain. It may be involved in regulation of cytoskeletal reorganization, cell proliferation, and cell motility, and acts as a tumor suppressor in hepatoma cells. The gene is located in a region of chromosome 13 that is associated with loss of heterozygosity in hepatocellular carcinomas. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Aug 2011]