

Product datasheet for **RC217832A1V**

Human NAA60 (NM_024845) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human NAA60 (NM_024845) AAV Particle
Tag:	Myc-DDK
Symbol:	NAA60
Synonyms:	HAT4; hNaa60; NAT15; NatF
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC217832 representing NM_024845 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGACAGAGGTGGTGCCATCCAGCGCGCTCAGCGAGGTCAGCCTGCGCCTCCTCTGCCACGATGACATAG
ACACTGTGAAGCACCTGTGTGGCGACTGGTTCCCCATCGAGTACCCAGACTCATGGATCGTGATATCAC
ATCCAACAAGAAGTCTTTTCCCTTGCTGCAACCTACAGAGGTGCCATTGTGGGAATGATAGTAGCTGAA
ATTAAGAACAGGACCAAAATACATAAAGAGGATGGAGATATTCTAGCATCCAACCTCTCTGTTGACACAC
AAGTCGCGTACATCCTAAGTCTGGGCGTCGTGAAAGAGTTCAGGAAGCACGGCATAGGTTCCCTCTTACT
TGAAAGTTTAAAGGATCACATATCAACCACCGCCAGGACCACTGCAAAGCCATTACCTGCATGTCCTC
ACCACCAACAACACAGCAATAAACTTCTATGAAAAACAGAGACTTCAAGCAGCACCACTATCTCCCTATT
ACTACTCCATTGAGGGGTCTCAAAGATGGCTTACCTATGTCCTCTACATCAACGGCGGGCACCCCTCC
CTGGACGATTTTGGACTACATCCAGCACCTGGGCTCTGCACTAGCCAGCCTGAGCCCCTGCTCCATTCCG
CACAGAGTCTACCGCCAGGCCACAGCCTGCTCTGCAGCTTCTGCCATGGTCGGGCATCTCTTCCAAGA
GTGGCATCGAGTACAGCCGACCATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC217832 representing NM_024845 Red=Cloning site Green=Tags(s) MTEVVPSSALSEVSLRLLCHDDIDTVKHLCDWFPPIEYPSWYRDITSNKKFFSLAATYRGAIVGMIVAE IKNRTKIHKEDGDILASNFSVDTQVAYILSLGVVKEFRKHGIGSLLLES�KDHISTTAQDHCKAIYLHVL TTNNTAINFYENRDFKQHHYLPYYYSIRGVLDGFTYVLYINGGHPWTILDYIQHLGSALASLSPCSIP HRVYRQAHSLLCSFLPWSGISSKSGIEYSRTM TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_024845
ORF Size:	726 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_024845.1, NP_079121.1</u>
RefSeq Size:	2553 bp
RefSeq ORF:	729 bp
Locus ID:	79903
UniProt ID:	<u>Q9H7X0</u>
Cytogenetics:	16p13.3
MW:	27.3 kDa