

Product datasheet for **RC217845L2V**

CLEC9A (NM_207345) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CLEC9A (NM_207345) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CLEC9A
Synonyms:	CD370; DNGR-1; DNGR1; UNQ9341
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_207345
ORF Size:	723 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217845).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_207345.1
RefSeq Size:	1674 bp
RefSeq ORF:	726 bp
Locus ID:	283420
UniProt ID:	Q6UXN8
Cytogenetics:	12p13.2
Protein Families:	Transmembrane
MW:	27.1 kDa



[View online »](#)

Gene Summary:

CLEC9A is a group V C-type lectin-like receptor (CLR) that functions as an activation receptor and is expressed on myeloid lineage cells (Huysamen et al., 2008 [PubMed 18408006]). [supplied by OMIM, Aug 2008]