

## Product datasheet for **RC217850L3V**

### POGZ (NM\_145796) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	POGZ (NM_145796) Human Tagged ORF Clone Lentiviral Particle
Symbol:	POGZ
Synonyms:	MRD37; WHSUS; ZNF280E; ZNF635; ZNF635m
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_145796
ORF Size:	3945 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217850).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_145796.3</a>
RefSeq Size:	6351 bp
RefSeq ORF:	3948 bp
Locus ID:	23126
UniProt ID:	<a href="#">Q7Z3K3</a>
Cytogenetics:	1q21.3
Domains:	zf-C2H2
MW:	145.3 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene appears to be a zinc finger protein containing a transposase domain at the C-terminus. This protein was found to interact with the transcription factor SP1 in a yeast two-hybrid system. Alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Aug 2010]