

Product datasheet for **RC217870L3V**

YIPF5 (NM_001024947) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	YIPF5 (NM_001024947) Human Tagged ORF Clone Lentiviral Particle
Symbol:	YIPF5
Synonyms:	FinGER5; SB140; SMAP-5; SMAP5; YIP1A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001024947
ORF Size:	771 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217870).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001024947.1
RefSeq Size:	3463 bp
RefSeq ORF:	774 bp
Locus ID:	81555
UniProt ID:	Q969M3
Cytogenetics:	5q31.3
Protein Families:	Transmembrane
MW:	28 kDa



[View online »](#)

Gene Summary: Plays a role in transport between endoplasmic reticulum and Golgi.[UniProtKB/Swiss-Prot Function]