

Product datasheet for **RC217875**

CD3D (NM_001040651) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: CD3D (NM_001040651) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: CD3D
Synonyms: CD3-DELTA; IMD19; T3D
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC217875 representing NM_001040651
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**C

ATGGAACATAGCAGTTTCTCTCTGGCCTGGTACTGGCTACCTTCTCTCGCAAGTGAGCCCTTCAAGA
TACCTATAGAGGAACCTGAGGACAGAGTGTGTGAATTGCAATACCAGCATCACATGGGTAGAGGGAAC
GGTGGGAACACTGCTCTCAGACATTACAAGACTGGACCTGGGAAAACGCATCCTGGACCCACGAGGAATA
TATAGGTGAATGGGACAGATATATAACAAGGACAAAGAATCTACCGTGCAAGTTCATTATCGAACTGCCG
ACACACAAGCTCTGTTGAGGAATGACCAGGTCTATCAGCCCTCCGAGATCGAGATGATGCTCAGTACAG
CCACCTTGGAGGAACTGGGCTCGGAACAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC217875 representing NM_001040651
Red=Cloning site **Green**=Tags(s)

MEHSTFLSGLVLATLLSQVSPFKIPIEELEDVFNVCNTSITWVEGTVGTLLSDITRLDLGKRILDPRGI
YRCNGTDIYKDKESTVQVHYRTADTQALLRNDQVYQPLRDRDDAQYSHLGGNWARNK

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001040651

ORF Size: 381 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

| | |
|-------------------------------|---|
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_001040651.2</u> |
| RefSeq Size: | 639 bp |
| RefSeq ORF: | 384 bp |
| Locus ID: | 915 |
| UniProt ID: | <u>P04234</u> |
| Cytogenetics: | 11q23.3 |
| Protein Families: | Druggable Genome |
| Protein Pathways: | Hematopoietic cell lineage, Primary immunodeficiency, T cell receptor signaling pathway |
| MW: | 14.48 kDa |
| Gene Summary: | <p>The protein encoded by this gene is part of the T-cell receptor/CD3 complex (TCR/CD3 complex) and is involved in T-cell development and signal transduction. The encoded membrane protein represents the delta subunit of the CD3 complex, and along with four other CD3 subunits, binds either TCR alpha/beta or TCR gamma/delta to form the TCR/CD3 complex on the surface of T-cells. Defects in this gene are a cause of severe combined immunodeficiency autosomal recessive T-cell-negative/B-cell-positive/NK-cell-positive (SCIDBNK). Two transcript variants encoding different isoforms have been found for this gene. Other variants may also exist, but the full-length natures of their transcripts has yet to be defined. [provided by RefSeq, Feb 2009]</p> |