

## Product datasheet for **RC217890L3V**

### EIF5AL1 (NM\_001099692) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	EIF5AL1 (NM_001099692) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EIF5AL1
Synonyms:	bA342M3.3; eIF-4D; eIF-5A; eIF-5A-1; eIF-5A1; EIF5A; EIF5AP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001099692
ORF Size:	462 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217890).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001099692.1</a>
RefSeq Size:	3849 bp
RefSeq ORF:	465 bp
Locus ID:	143244
UniProt ID:	<a href="#">Q6IS14</a>
Cytogenetics:	10q22.3
MW:	16.6 kDa



[View online »](#)

**Gene Summary:**

mRNA-binding protein involved in translation elongation. Has an important function at the level of mRNA turnover, probably acting downstream of decapping. Involved in actin dynamics and cell cycle progression, mRNA decay and probably in a pathway involved in stress response and maintenance of cell wall integrity. Functions as a regulator of apoptosis. Mediates effects of polyamines on neuronal process extension and survival. May play an important role in brain development and function, and in skeletal muscle stem cell differentiation (By similarity).[UniProtKB/Swiss-Prot Function]