

Product datasheet for **RC217893**

C17orf77 (NM_152460) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C17orf77 (NM_152460) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	C17orf77
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC217893 representing NM_152460 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGATGAATTGGCGCTGTCATTTCCCTGACATGTCTCCTTCCTGAGAACAGAGCGTCTGTCTCCAA
 GCCAGCCGCTGTCATTCCAGTGCCTTAAAGCTCCAGCCCACTCACCTGGGAGGATGAGAAGCAGCAGAG
 ATGGGGACAGCCACATGGTCCTGTCTCCTCTCCTCTCCTTGGTGATCACAGGTGTCTGGTGCCCTTCAGG
 GACTTGAATCCAAGCTCGGAAGTCAACACCGCAAACCTACTCGAATCTCCTTCAAGTCTCTTGCTCACTT
 CCTGTTACATCTGTAGTACTTCTCTTTTACATATTAGGAGAAAAGAGGTGTCATTCTTTGAAAAGACT
 AAGGTACAGTGTATGTTGCAAGGTGTGCCAAACTTTTGTGCATGTGAAAAGAAAATGTGTCCGGGACT
 GGTCAAGTTTGCACGGGAGTGCACGTAGGAGCCAAGGAGCAGGAGGCCAGGAGGAACCAAGCATTGA
 GGAGCTGCGGGATCTACTGCCTGGAGGAGAGGACGGACAAGGCCAGCCATGAGGAGTGCAGGGAGAGGAG
 CACTTTGGGAAGACCACAGTGCACAGGACTCACACCCAGCCTTGCTGGAGAAAGCCCATGCCCCAGGCTT
 CTTCCAGGCAGCCCCACTGTAAGACATCTATTGCCTCTTCTGCCCTGGGCTGTCTGACCCACTGCCCC
 TTCCCCAGGAACACTTCCACTGGGATCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >RC217893 representing NM_152460
 Red=Cloning site Green=Tags(s)


MDELALSFSLTCLLPENRASLSPSQPLSFQCLKAPATLTWEDEKQQRWGQPHGPVSSPLLGDHRCCLVPFR
 DLNPSSEVNTANLLESPSSLLL TSCYICSYFSFYILGEKRCHSLKRLRYSVCCCKVCPNFCACGKENVSGT
 GQVCTGVHVGAKQEPEGGTQALRSCGIYCLEERTDKASHEECRERSTLGRPQCTGLTPSLAGESPCPRL
 LPGSPTVRHLIASSCPGLSDPLPLPPGTLPLGS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8012_e05.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_152460

ORF Size: 729 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_152460.2](#), [NP_689673.2](#)

RefSeq Size: 2525 bp

RefSeq ORF: 732 bp

Locus ID: 146723


UniProt ID: [Q96MU5](#)

Cytogenetics: 17q25.1


Protein Families: Secreted Protein

MW: 26.1 kDa

Product images:



Circular map for RC217893



Western blot validation of overexpression lysate (Cat# [LY407520]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC217893 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).