

Product datasheet for **RC217895L3V**

MSL3L1 (MSL3) (NM_078628) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MSL3L1 (MSL3) (NM_078628) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MSL3L1
Synonyms:	MRSXBA; MRXS36; MRXSBA; MSL3L1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_078628
ORF Size:	1248 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217895).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_078628.1 , NP_523352.1
RefSeq Size:	3564 bp
RefSeq ORF:	1251 bp
Locus ID:	10943
UniProt ID:	Q8N5Y2
Cytogenetics:	Xp22.2
Domains:	CHROMO
Protein Families:	Transcription Factors



[View online »](#)

MW: 47.8 kDa

Gene Summary: This gene encodes a nuclear protein that is similar to the product of the Drosophila male-specific lethal-3 gene. The Drosophila protein plays a critical role in a dosage-compensation pathway, which equalizes X-linked gene expression in males and females. Thus, the human protein is thought to play a similar function in chromatin remodeling and transcriptional regulation, and it has been found as part of a complex that is responsible for histone H4 lysine-16 acetylation. This gene can undergo X inactivation. Alternative splicing results in multiple transcript variants. Related pseudogenes have been identified on chromosomes 2, 7 and 8. [provided by RefSeq, Jul 2010]