

Product datasheet for RC217906

OR10Z1 (NM_001004478) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OR10Z1 (NM_001004478) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	OR10Z1
Synonyms:	OR1-15
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC217906 representing NM_001004478 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGGCAGACCAACGTAACCTCCTGGAGGGATTTGTCTTCCTGGGCTTCTCCAGTTCTGGGGAGTTGC
AGCTCCTTCTTTGCCTTGTTCCTCTCTGTATCTAGTCACTCTGACCAGCAATGTCTTCATTATCAT
AGCCATCAGGCTGGATAGCCATCTGCACACCCCATGTACCTCTCCTTTCCTTCTATCCTTCTCTGAG
ACCTGCTACACTTTGGGCATCATCCCTAGAATGCTCTCTGGCCTGGCTGGGGGGACCAGGCTATCTCT
ATGTGGGCTGTGCTGCCAGATGTTCTTTCTGCCTCATGGCCTGTAATACTGCTTCTTCTGGCTGC
CATGGGCTTTGACAGATATGTGGCCATCTGTGCTCCACTCCACTATGCCAGCCACATGAATCCTACCCTC
TGTGCCAGCTGGTCATTACTTCTTCTGACTGGATACCTCTTTGGACTGGGAATGACACTAGTATTT
TCCACCTCTATTCTGCAGCTCCCATGAAATCCAGCACTTTTTTTGTGACACGCCACCTGTGCTGAGCCT
AGCCTGTGGAGATACAGGCCGAGTGAGCTGAGGATCTTATCCTCAGTCTTTGGTCTCTTGGTCTCC
TTCTTCTCATCACCATCTCTACGCCTACATCTTGGCAGCAATACTGAGGATCCCCTCTGCTGAGGGGC
AGAAGAAGGCCTTCTCCACTTGTGCTCGCACCTTACAGTGGTCATTATTCATTATGGCTGTGCTTCTT
CGTGTACCTGAGGCCAAAGCCAGCTACTCTCTGAGAGAGATCAGCTTATTGCCATGACCTATACTGTA
GTGACCCCTCCTTAATCCATTGTTTATAGTCTAAGGAATAGGGCTATACAGACAGCTCTGAGGAATG
CTTTCAGAGGGAGATTGCTGGTAAAGGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC217906 representing NM_001004478
Red=Cloning site Green=Tags(s)

MGQTNVTSWRDFVFLGFSSSGELQLLLFALFLSLYLVTLSNVFIIIAIRLDShLHTPMYLFSLFSFSE
 TCYTLGIIPRML SGLAGGDQAI SYVGCAQMFFSASWACTNCFLLAAMGFDRYVAICAPLHYASHMNP TL
 CAQLVITSFLTGYL FGLGMLVIFHL SFCSSHEIQHFFCDTPPVLSLACGDTGPSELRI FLSLLVLLVS
 FFFITISYAYILAAILRIPSAEQKKAFSTCASHLTVV I IHYGCASFVYL RPKASYSLERDQLIAMTYTV
 VTPLLNP I VYSLRNRAIQ T ALRNAFRGRLLGKG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8002_c01.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001004478

ORF Size: 939 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001004478.1, NP_001004478.1</u>
RefSeq Size:	942 bp
RefSeq ORF:	942 bp
Locus ID:	128368
UniProt ID:	<u>Q8NGY1</u>
Cytogenetics:	1q23.1
Protein Families:	Transmembrane
Protein Pathways:	Olfactory transduction
MW:	34.4 kDa
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC217906



Western blot validation of overexpression lysate (Cat# [LY423781]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC217906 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).