

## Product datasheet for **RC217911L4V**

### SEC23B (NM\_006363) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SEC23B (NM_006363) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SEC23B
Synonyms:	CDA-II; CDAII; CDAN2; CWS7; HEMPAS; hSec23B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006363
ORF Size:	2301 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217911).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006363.4</a>
RefSeq Size:	3429 bp
RefSeq ORF:	2304 bp
Locus ID:	10483
UniProt ID:	<a href="#">Q15437</a>
Cytogenetics:	20p11.23
Domains:	zf-Sec23_Sec24, Sec23_trunk, Sec23_helical, Gelsolin
MW:	86.5 kDa


[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a member of the SEC23 subfamily of the SEC23/SEC24 family, which is involved in vesicle trafficking. The encoded protein has similarity to yeast Sec23p component of COPII. COPII is the coat protein complex responsible for vesicle budding from the ER. The function of this gene product has been implicated in cargo selection and concentration. Multiple alternatively spliced transcript variants have been identified in this gene. [provided by RefSeq, Feb 2010]