

## Product datasheet for **RC217912L4V**

### NT5C1B (NM\_001002006) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	NT5C1B (NM_001002006) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NT5C1B
Synonyms:	AIRP; CN-IB; CN1B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001002006
ORF Size:	1830 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217912).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001002006.1</a>
RefSeq Size:	2864 bp
RefSeq ORF:	1833 bp
Locus ID:	93034
UniProt ID:	<a href="#">Q96P26</a>
Cytogenetics:	2p24.2
Protein Pathways:	Metabolic pathways, Nicotinate and nicotinamide metabolism, Purine metabolism, Pyrimidine metabolism



[View online »](#)

**MW:** 68.6 kDa

**Gene Summary:** Cytosolic 5-prime nucleotidases, such as NT5C1B, catalyze production of adenosine, which regulates diverse physiologic processes (Sala-Newby and Newby, 2001 [PubMed 11690631]). [supplied by OMIM, Mar 2008]