

Product datasheet for **RC217919L4V**

AMPD1 (NM_000036) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	AMPD1 (NM_000036) Human Tagged ORF Clone Lentiviral Particle
Symbol:	AMPD1
Synonyms:	MAD; MADA; MMDD
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_000036
ORF Size:	2241 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217919).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000036.1 , NP_000027.1
RefSeq Size:	2426 bp
RefSeq ORF:	2244 bp
Locus ID:	270
UniProt ID:	P23109
Cytogenetics:	1p13.2
Domains:	A_deaminase
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Metabolic pathways, Purine metabolism

MW: 86.5 kDa

Gene Summary: Adenosine monophosphate deaminase 1 catalyzes the deamination of AMP to IMP in skeletal muscle and plays an important role in the purine nucleotide cycle. Two other genes have been identified, AMPD2 and AMPD3, for the liver- and erythrocyte-specific isoforms, respectively. Deficiency of the muscle-specific enzyme is apparently a common cause of exercise-induced myopathy and probably the most common cause of metabolic myopathy in the human. Alternatively spliced transcript variants encoding different isoforms have been identified in this gene.[provided by RefSeq, Feb 2010]