

Product datasheet for **RC217929A1V**

Human PLAUR (NM_001005376) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human PLAUR (NM_001005376) AAV Particle
Tag:	Myc-DDK
Symbol:	PLAUR
Synonyms:	CD87; U-PAR; UPAR; URKR
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC217929 representing NM_001005376 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGTCACCCGCCGCTGCTGCCGCTGCTGCTGCTGCTCCACACCTGCGTCCCAGCCTCTTGGGGCCTGC
GGTGCATGCAGTGTAAAGACCAACGGGGATTGCCGTGTGGAAGAGTGCGCCCTGGGACAGGACCTCTGCAG
GACCACGATCGTGCCTTGTGGGAAGAAGGAGAAGAGCTGGAGCTGGTGGAGAAAAGCTGTACCCACTCA
GAGAAGACCAACAGGACCCTGAGCTATCGGACTGGCTTGAAGATCACCAGCCTTACCGAGGTTGTGTGTG
GGTTAGACTTGTGCAACCAGGGCAACTCTGGCCGGGCTGTACCTATTCGGAAGCCGTTACCTCGAATG
CATTTCTGTGGCTCATCAGACATGAGCTGTGAGAGGGGCCGGCACCAGAGCCTGCAGTGCCGCAGCCCT
GAAGAACAGTGCCTGGATGTGGTGACCCACTGGATCCAGGAAGGTGAAGAAGGGCGTCCAAGGATGACC
GCCACCTCCGTGGCTGTGGCTACCTTCCCGGCTGCCCGGGCTCCAATGGTTTCCACAACAACGACACCTT
CCACTTCTGAAATGCTGCAACACCACCAAAATGCAACGAGGGCCCAATCCTGGAGCTTGAAAATCTGCCG
CAGAATGGCCGCCAGTGTACAGCTGCAAGGGGAACAGCACCCATGGATGCTCCTCTGAAGAGACTTTCC
TCATTGACTGCCGAGGCCCATGAATCAATGTCTGGTAGCCACCGGCACTACGAACGCTCACTCTGGGG
AAGCTGGTTGCCATGTAAAGTACTACTGCCCTGAGACCACCATGCTGTGAGGAAGCCCAAGCTACTCAT
GTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC217929 representing NM_001005376 Red=Cloning site Green=Tags(s) MGHPPLLPLLLLHTCVPASWGLRCMQCKTNGDCRVEECALGQDLCRTTIVRLWEEGEELELVEKSCTHS EKTNRTLSYRTGLKITSLTEVVCGLDLNQGNSGRAVTYSRSRYLECISGSSDMSCERGRHQLQCRSP EEQCLDVVTHWIQEGEEGRPKDDRHLRGCGYLPGCPGSGNFHNNDTFHFLKCCNTTKCNEGPILELENLP QNGRQCYSCKGNSTHGCSSEETFLIDCRGPMNQCLVATGOTHERSLWGSWLPCKSTTALRPPCCEEAQATH V TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001005376
ORF Size:	843 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001005376.1</u>
RefSeq Size:	1437 bp
RefSeq ORF:	846 bp
Locus ID:	5329
UniProt ID:	<u>Q03405</u>
Cytogenetics:	19q13.31
MW:	31.26 kDa