

Product datasheet for **RC217941L2V**

CD200 (NM_001004196) Human Tagged ORF Clone Lentiviral Particle

Product data:


Product Type:	Lentiviral Particles
Product Name:	CD200 (NM_001004196) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CD200
Synonyms:	MOX1; MOX2; MRC; OX-2
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001004196
ORF Size:	882 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217941).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001004196.1 , NP_001004196.1
RefSeq Size:	2247 bp
RefSeq ORF:	885 bp
Locus ID:	4345
UniProt ID:	P41217
Cytogenetics:	3q13.2
Protein Families:	Transmembrane
MW:	32.8 kDa




[View online »](#)

Gene Summary:

This gene encodes a type I membrane glycoprotein containing two extracellular immunoglobulin domains, a transmembrane and a cytoplasmic domain. This gene is expressed by various cell types, including B cells, a subset of T cells, thymocytes, endothelial cells, and neurons. The encoded protein plays an important role in immunosuppression and regulation of anti-tumor activity. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jan 2016]

Product images:

[RC217941L2] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC217941L2V particle to overexpress human CD200-mGFP fusion protein.



[RC217941L2] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC217941L2V particle to overexpress human CD200-mGFP fusion protein.