

## Product datasheet for RC217961L3V

### CTNNAL1 (NM\_003798) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CTNNAL1 (NM_003798) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CTNNAL1
Synonyms:	ACRP; alpha-CATU; CLLP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_003798
ORF Size:	2202 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217961).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_003798.1</a>
RefSeq Size:	2573 bp
RefSeq ORF:	2205 bp
Locus ID:	8727
UniProt ID:	<a href="#">Q9UBT7</a>
Cytogenetics:	9q31.3
Protein Families:	Druggable Genome
MW:	81.9 kDa


[View online »](#)

**Gene Summary:**

May modulate the Rho pathway signaling by providing a scaffold for the Lbc Rho guanine nucleotide exchange factor (ARHGEF1).[UniProtKB/Swiss-Prot Function]