

Product datasheet for RC217981

Arp2 (ACTR2) (NM_001005386) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Arp2 (ACTR2) (NM_001005386) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Arp2
Synonyms:	ARP2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC217981 representing NM_001005386 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACAGCCAGGGCAGGAAGGTGGTGGTGTGCGACAACGGCACCGGTTTGTGAAGTGTGGATATGCAG
GCTCTAACTTTCCAGAACACATCTCCAGCTTTGGTTGGAAGACCTATTATCAGATCAACCACCAAAGT
GGGAAACATTGAAATCAAGAATAACAAAAAGATGGATCTTATGGTTGGTGTGATGAGGCAAGTGAATTACGA
TCAATGTTAGAAGTTAACTACCCTATGGAAAATGGCATAGTACGAAATTGGGATGACATGAAACACCTGT
GGGACTACACATTTGGACCAGAGAACTTAATATAGATACCAGAAATTGAAAATCTTACTCACAGAACC
TCCTATGAACCAACCAAAAACAGAGAGAAGATTGTAGAGGTAATGTTTGAAACTTACCAGTTTTCCGGT
GTATATGTAGCCATCCAGGCAGTTCTGACTTTGTACGCTCAAGGTTTATTGACTGGTGTAGTGGTAGACT
CTGGAGATGGTGTGACTCACATTTGCCAGTATATGAAGGCTTTCTCTCCCTCATCTTACCAGGAGACT
GGATATTGCTGGGAGGGATATAACTAGATATCTTATCAAGCTACTTCTGTTGCGAGGATACGCCTTCAAC
CACTCTGCTGATTTTGAACGGTTCGCATGATTAAGAAAAACTGTGTTACGTGGGATATAATATTGAGC
AAGAGCAGAACTGGCCTTAGAAACCACAGTATTAGTTGAATCTTATACACTCCCAGATGGACGTATCAT
CAAAGTTGGGGGAGAGAGATTTGAAGCACCAGAAGCTTTATTTACGCCTCACTTGATCAATGTTGAAGGA
GTTGGTGTGCTGAATTGCTTTTAAACACAATTCAGGCAGTGTGATACCAGATCTGAATTTACA
AACACATTGTGCTTTCTGGAGGGTCTACTATGTATCCTGCCTGCATCACGGTTGGAACGAGAACTTAA
ACAGCTTTACTTAGAACGAGTTTTGAAGGGTGTGTTGAAAAACTTTCTAAATTTAAGATCCGCATTGAA
GACCCACCCCGCAGAAAGCACATGGTATTCCTGGGTGGTGCAGTTCTAGCGGATATCATGAAAGACAAAG
ACAACCTTTGGATGACCCGACAAGAGTACCAAGAAAAGGGTGTCCGTGTGCTAGAGAACTTGGTGTGAC
TGTTCGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RC217981 representing NM_001005386
Red=Cloning site Green=Tags(s)

MDSQGRKVVVCDNGTGFVKCGYAGSNFPEHIFPALVGRPIIRSTTKVGNIEIKNNKKMDLMVGDEASELR
 SMLEVNYPMENGIVRNWDDMKHLWDYTFGPEKLNIDTRNCKILLTEPPMNPTKNREKIVEVMFETYQFSG
 VYVAIQAVLTLYAQGLLTGVVVDSGDGVTHICPVYEGFSLPHLTRLLDIAGRDITRYLIKLLLRGYAFN
 HSADFETVRMIKEKLCYVGYNIEQEQLALETTVLVESYTLPDGRIKIVGGERFEAPEALFQPHLINVEG
 VGVAELLFNTIQAADIDTRSEFYKHIVLGGSTMYPLPSRLERELKQLYLERVLKGDVEKLSKFKIRIE
 DPPRRKHMVFLGGAVLADIMKDKDNFWMTRQEYQEKGVRVLEKLGVTVR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6498_c11.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001005386

ORF Size: 1197 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001005386.2](#)

RefSeq Size: 3878 bp

RefSeq ORF: 1200 bp

Locus ID: 10097

UniProt ID: [P61160](#)


Cytogenetics: 2p14

Protein Families: Druggable Genome


MW: 45.2 kDa

Gene Summary: The specific function of this gene has not yet been determined; however, the protein it encodes is known to be a major constituent of the ARP2/3 complex. This complex is located at the cell surface and is essential to cell shape and motility through lamellipodial actin assembly and protrusion. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC217981



Western blot validation of overexpression lysate (Cat# [LY423711]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC217981 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).