

Product datasheet for **RC217995A1V**

Human hnRNP F (HNRNPF) (NM_001098206) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human hnRNP F (HNRNPF) (NM_001098206) AAV Particle
Tag:	Myc-DDK
Symbol:	hnRNP F
Synonyms:	HNRPF; mcs94-1; OK/SW-cl.23
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC217995 representing NM_001098206 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGATGCTGGGCCCTGAGGGAGGTGAAGGCTTTGTGGTCAAGCTCCGTGGCCTGCCCTGGTCTGCTCTG
TTGAGGACGTGCAGAATTCTCTCTGACTGCACGATTCATGATGGGGCCGAGGTGTCCATTTTCATCTA
CACTAGAGAGGGCAGGCAGAGTGGTGAAGCTTTTGTGAACCTGGATCAGAAGATGATGTAAAAATGGCC
CTGAAAAAAGACAGGGAAGCATGGGACACCGGTACATTGAGGTGTTGAGTCCCACAGAACCGAGATGG
ATTGGGTGTTGAAGCACAGTGGTCCCAACAGTGCCGACAGCGCCAACGATGGCTTCGTGCGGCTTCGAGG
ACTCCCATTTGGATGCACAAAGGAAGAAATTGTTCAAGTTCTTCTCAGGGTTGGAATTTGTGCCAAACGGG
ATCACATTGCCTGTGGACCCCGAAGGCAAGATTACAGGGGAAGCGTTCGTGCAGTTTGCCTCGCAGGAGT
TAGCTGAGAAGGCTCTAGGGAACACAAGGAGAGGATAGGGCACAGGTACATTGAGGTGTTAAGAGCAG
CCAGGAGGAAGTTAGGTCATACTCAGATCCCCCTCTGAAGTTCATGTCCGTGCAGCGGCCAGGGCCCTAT
GACCGGCCCGGGACTGCCAGGAGGTACATTGGCATCGTGAAGCAGGCAGGCCTGGAAAGGATGAGGCCTG
GTGCCTACAGCACAGGCTACGGGGGCTACGAGGAGTACAGTGGCCTCAGTGATGGCTACGGCTTCACCAC
CGACCTGTTCCGGAGAGACCTCAGCTACTGTCTCTCCGGAATGTATGACCACAGATACGGCGACAGTGAG
TTCACAGTGCAGAGCACACAGGCCACTGTGTCCACATGAGGGGCCTGCCGTACAAAGCGACCGAGAACG
ACATTTACAATTCTTCTCTCTCTCAACCCTGTGAGAGTCCATATTGAGATTGGCCAGATGGAAGAGT
GACGGGTGAAGCAGATGTTGAGTTTGCTACTCATGAAGAAGCTGTGGCAGCTATGTCAAAGACAGGGCC
AATATGCAGCACAGATATATAGAACTCTTCTTGAATTCACAACAGGGGCCAGCAATGGGGCGTATAGCA
GCCAGGTGATGCAAGGCATGGGGGTGTCTGCTGCCAGGCCACTTACAGTGGCCTGGAGAGCCAGTCAGT
GAGTGGCTGTTACGGGGCCGGCTACAGTGGGCAGAACAGCATGGGTGGCTATGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC217995 representing NM_001098206 Red=Cloning site Green=Tags(s) MMLGPEGGEGFVVKLRGLPWSCSVEDVQNFLSDCTIHDGAAGVHFIYTREGRQSGEAFVELGSEDDVKMA LKKDRESMGHRYIEVFRSHRTEMDWVLKHSGPNSADSANDGFVRLRGLPFGCTKEEIVQFFSGLEIVPNG ITLPVDPGKITGEAFVQFASQELA EKALGKHKERIGHRYIEVFKSSQEEVRSYSDPPLKFMSVQRP GPY DRPGTARRYIGIVKQAGLERMRPGAYSTGYGGYEEYSGLS DGYGFTTDLFGRDLSYCLSGMYDHRYGDSE FTVQSTTGHCVMRGLPYKATENDIYNFFSPLNPVRVHIEIGPDGRVTGEADVEFATHEEAVAAMSKDRA NMQHRYIELFLNSTTGASNGAYSSQVMQGMGVSAQA TYSGLESQSVSGCYGAGYSGQNSMGGYD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001098206
ORF Size:	1245 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001098206.1</u> , <u>NP_001091676.1</u>
RefSeq Size:	2667 bp
RefSeq ORF:	1248 bp
Locus ID:	3185
UniProt ID:	<u>P52597</u>
Cytogenetics:	10q11.21
MW:	45.7 kDa