

Product datasheet for **RC218002L3V**

CDX4 (NM_005193) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CDX4 (NM_005193) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CDX4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_005193
ORF Size:	852 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218002).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005193.1 , NP_005184.1
RefSeq Size:	855 bp
RefSeq ORF:	855 bp
Locus ID:	1046
UniProt ID:	O14627
Cytogenetics:	Xq13.2
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transcription Factors
MW:	30.3 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of a small subfamily of homeobox containing transcription factors. The encoded protein may regulate homeobox gene expression during anteroposterior patterning and hematopoiesis. [provided by RefSeq, Aug 2012]