

Product datasheet for **RC218064**

LILRB3 (NM_001081450) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LILRB3 (NM_001081450) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	LILRB3
Synonyms:	CD85A; HL9; ILT-5; ILT5; LILRA6; LIR-3; LIR3; PIR-B; PIRB
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC218064 representing NM_001081450
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGACGCCCGCCCTCACAGCCCTGCTCTGCCTTGGGCTGAGTCTGGGCCCCAGGACCCCGTGCAGGCAG
GGCCCTTCCCAAACCCACCCCTCTGGGCTGAGCCAGGCTCTGTGATCAGCTGGGGAGCCCGTGACCAT
CTGGTGTGAGGGGAGCCAGGAGGCCAGGAGTACCAGACTGCATAAAGAGGGAAGCCAGAGCCCTTGGAC
AGAAATAACCCACTGGAACCCAAGAACAAGGCCAGATTCTCCATCCCATCCATGACAGAGCACCATGCAG
GGAGATACCGCTGCCACTATTACAGCTCTGAGGCTGGTCCAGAGCCAGCGACCCCTGGAGATGGTGTGAT
GACAGGAGCTACAGCAAACCCACCCCTCTCAGCCCTGCCAGCCCTGTGGTGGCCCTCAGGGGGAATATG
ACCCTCCGATGTGGCTCACAGAAGGGATATCACCATTTTGTCTGATGAAGGAAGGAGAACCAGCTCC
CCCGGACCCTGGACTCACAGCAGCTCCACAGTCGGGGTCCAGGCCCTGTTCCCTGTGGCCCCGTGAC
CCCCAGCCACAGGTGGAGTTTACATGCTATTACTATTATACAAACCCCTGGGTGTGGTCCCACCCC
AGTGACCCCTGGAGATTCTGCCCTCAGGCGTGTCTAGGAAGCCCTCCCTCCTGACCCTGCAGGGCCCTG
TCCTGGCCCTGGGCAGAGCCTGACCCTCCAGTGTGGCTCTGATGTGGGTACAACAGATTTGTTCTGTA
TAAGGAGGGGAAACGTGACTTCCCTCCAGCGCCCTGGCCAGCAGCCCAGGCTGGGCTCTCCAGGCCAAC
TTCACCCTGGGCCCTGTGAGCCCTCCAATGGGGCCAGTACAGGTGCTACGGTGACACAACTCTCCT
CCGAGTGGTGGCCCCCAGCGACCCCTGAACATCCTGATGGCAGGACAGATCTATGACACCGTCTCCCT
GTCAGCACAGCCGGGCCACAGTGGCCTCAGGAGAGAACGTGACCCTGCTGTGTGATCATGGTGGCAG
TTTGACACTTTCCTTCTGACCAAAGAAGGGCAGCCATCCCCACTGCGTCTGAGATCAATGTACGGAG
CTCATAAGTACCAGGCTGAATTCCCATGAGTCTGTGACCTCAGCCACGCGGGGACCTACAGGTGTA
CGGCTCATACAGCTCAAACCCACCTGCTGTCTCACCCAGTGAGCCCTGGAGCTCGTGGTCTCAGGA
CACTCTGGAGGCTCCAGCCTCCCACCCACAGGGCCGCTCCACACCTGGTCTGGGAAGATACCTGGAGG
TTTTGATTGGGGTCTCGGTGGCCTTCGTCTGCTCTTCTCCTCCTCTTCTCCTCCTCCGACGTCA
GCGTCACAGCAAACACAGGACATCTGACCAGAGAAAGACTGATTTCCAGCGTCTGCAGGGGCTGCGGAG
ACAGAGCCCAAGGACAGGGCCTGCTGAGGAGGTCCAGCCAGCTGCTGACGTCCAGGAAGAAAACCTCT
ATGCTGCCGTGAAGGACACAGTCTGAGGACAGGTTGGAGCTGGACAGTCCAGAGCCACACGATGA
AGACCCCAAGGACAGTATGCCCGGTGAAACACTCCAGTCTAGGAGAGAAATGGCCTCTCCTCCC
TCCTACTGTCTGGGAATTCCTGGACACAAGGACAGACAGGTGGAAGAGGACAGGAGATGGACTG
AGGCTGCTGCATCTGAAGCCTCCAGGATGTGACCTACGCCAGCTGCACAGCTTACCCTTAGACGGAA
GGCAACTGAGCCTCCTCCATCCCAGGAAGGGAACTCCAGCTGAGCCAGCATCTACGCCACTCTGGCC
ATCCAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTAA

Protein Sequence: >RC218064 representing NM_001081450
Red=Cloning site Green=Tags(s)

```
MTPALTALLCLGLSLGPRTRVQAGFPFKPTLWAEPGSVISWGSPVTIWCQGSQEAQEYRLHKEGSPEPLD
RNNPLEPKNKARFSIPSMTEHHAGRYRCHYYSSAGWSEPSDPLEMVMGTAYSKPTLSALPSVVASGGNM
TLRCGSQKGYHHFVLMKEGEHQLPRTLDSQQLHSRGFQALFPVGPVTPSHRWRF T C Y Y Y Y T N T P W V W S H P
SDPLEIILPSGVSRKPSLLTLQGPVLAPGQSLTLQCGSDVGYNRFVLYKEGERDFLQRPQQPQAGLSQAN
FTLGPVSPSNGGQYRCYGAHNLSSSEWAPSDDLNILMAGQIYDTVSLSAQPGPTVASGENVTLLCQSWWQ
FDTFLLTKEGAAHPPLRLRSMYGAKHYQAEFPMSPVTSAHAGTYRCYGSYSSNPHELLSHPSEPELELVVSG
HSGGSSLPTGPPSTPGLGRYLEVLIGVSAFVLLLLFLLLFLRRQRHSHKRTSDQRKTD F Q R P A G A A E
TEPKDRGLLRSSPAADVQEENLYAAVKDTQSEDRVELDSQQSPHDEDPQAVTYAPVKHSSPREMASPP
SSLSGEFLDKDRQVEEDRQMDTEAAASEASQDVTYAQLHSLTLRRKATEPPPSQEGEPPAEP SIYATLA
IH
```

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8070_f04.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001081450

ORF Size: 1896 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001081450.1, NP_001074919.1</u>
RefSeq Size:	2843 bp
RefSeq ORF:	1899 bp
Locus ID:	11025
Cytogenetics:	19q13.42
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	B cell receptor signaling pathway
MW:	69.3 kDa
Gene Summary:	<p>This gene is a member of the leukocyte immunoglobulin-like receptor (LIR) family, which is found in a gene cluster at chromosomal region 19q13.4. The encoded protein belongs to the subfamily B class of LIR receptors which contain two or four extracellular immunoglobulin domains, a transmembrane domain, and two to four cytoplasmic immunoreceptor tyrosine-based inhibitory motifs (ITIMs). The receptor is expressed on immune cells where it binds to MHC class I molecules on antigen-presenting cells and transduces a negative signal that inhibits stimulation of an immune response. It is thought to control inflammatory responses and cytotoxicity to help focus the immune response and limit autoreactivity. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p>

Product images:



Circular map for RC218064



Western blot validation of overexpression lysate (Cat# [LY421154]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC218064 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).