

Product datasheet for **RC218080A1V**

Human FGFR1OP2 (NM_015633) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human FGFR1OP2 (NM_015633) AAV Particle
Tag:	Myc-DDK
Symbol:	FGFR1OP2
Synonyms:	HSPC123-like; WIT3.0
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC218080 representing NM_015633 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAGTTGCACAATTGAGAAGGCACTTGCCGACGCTAAAGCTCTTGTTGAAAGATTAAGAGATCATGACG
ATGCAGCAGAATCTCTGATTGAGCAAACACAGCTCTCAACAAGCGAGTAGAAGCCATGAAACAGTATCA
GGAAGAAATTCAGAACCTAATGAAGTCGCGAGACATCGGCCACGGTCCACGTTAGTTATGGGAATCCAG
CAAGAAAACAGACAAATCAGAGAGTTGCAACAAGAAAACAAAGAATTACGTACATCTCTGGAAGAATC
AGTCGGCCTTGGAACCTATAATGAGCAAGTACCGAGAACAATGTTTAGATTGCTAATGGCTAGCAAAAA
AGATGATCCGGGTATAATAATGAAGTTAAAAGAGCAGCACTCCAAGATTGACATGGTACATCGTAACAAG
TCCGAAGGATTCTTCTTGATGCATCTCGACACATCCTTGAAGCACCTCAACATGGACTGGAGAGAAGGC
ACTTGGAAGCAAATCAGAAATGAAGTGAAGCACATGTTGACCAGATAACTGAAATGGCAGCAGTAATGAG
GAAAGCCATTGAAATTGACGAGCAACAGGGTTGCAAGGAACAAGAACGAATATTTCAACTGAACAAGAA
AACAAAGGCTTGAGAGAGATCCTTCAAATAACTCGAGAATCATTTTTGAACCTAAGGAAAGATGATGCGT
CGGAAAGTACTTCTTTGTCAGCATTAGTGACCAACAGTGATTTGAGTCTGAGGAAGAGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC218080 representing NM_015633 Red=Cloning site Green=Tags(s) MSCTIEKALADAKALVERLRDHDAAESLIEQTTALNKRVEAMKQYQEEIQELNEVARHRPRSTLVMGIQ QENRQIRELQQENKELRTSLEEHQSALELIMSKYREQMFRLLMASKKDDPGIIMKLKEQHSKIDMVHRNK SEGFFLDASRHILEAPQHGLERRHLEANQNELQAHVDQITEMAAYMRKAIEIDEQQGCKEQERIFQLEQE NKGLREILQITRESFLNLRKDDASESTLSALVTNSDLSLRKS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_015633
ORF Size:	759 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_015633.1</u>
RefSeq Size:	2306 bp
RefSeq ORF:	762 bp
Locus ID:	26127
UniProt ID:	<u>Q9NVK5</u>
Cytogenetics:	12p11.23
MW:	29.2 kDa