

Product datasheet for **RC218089L3V**

TTC8 (NM_144596) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TTC8 (NM_144596) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TTC8
Synonyms:	BBS8; RP51
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_144596
ORF Size:	1545 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218089).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_144596.2
RefSeq Size:	2248 bp
RefSeq ORF:	1548 bp
Locus ID:	123016
UniProt ID:	Q8TAM2
Cytogenetics:	14q31.3
Domains:	TPR
MW:	58.2 kDa



[View online »](#)

Gene Summary:

This gene encodes a protein that has been directly linked to Bardet-Biedl syndrome. The primary features of this syndrome include retinal dystrophy, obesity, polydactyly, renal abnormalities and learning disabilities. Experimentation in non-human eukaryotes suggests that this gene is expressed in ciliated cells and that it is involved in the formation of cilia. A mutation in this gene has also been implicated in nonsyndromic retinitis pigmentosa. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]