

## Product datasheet for RC218114A1V

### Human OR4D11 (NM\_001004706) AAV Particle

#### Product data:

|                           |                                                                                |
|---------------------------|--------------------------------------------------------------------------------|
| Product Type:             | AAV Particles                                                                  |
| Product Name:             | Human OR4D11 (NM_001004706) AAV Particle                                       |
| Tag:                      | Myc-DDK                                                                        |
| Symbol:                   | OR4D11                                                                         |
| Synonyms:                 | OR4D11P                                                                        |
| Mammalian Cell Selection: | None                                                                           |
| Vector:                   | pAAV-AC-Myc-DDK (PS100089)                                                     |
| ORF Nucleotide Sequence:  | >RC218114 representing NM_001004706<br>Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**

ATGGAGTTGGGAAATGTCACCAGAGTAAAAGAATTTATTTCTGGGACTTACTCAATCCCAAGACCAGA  
 GTTTGGTCTTGTCTTTTTTATGTCTTGTTACATGACGACTCTGCTGGGAAACCTCCTCATCATGGT  
 CACCGTGACCTGTGAGTCTCGCCTTCACACCCCATGTACTTCCTGCTCCGCAATCTAGCCATCCTTGAC  
 ATCTGCTTCTCCTCCACAAGTCTCTAAAGTCTTGCTGGACCTTCTGTCAAAGAAAAAGACCATATCCT  
 ATACAAGCTGCATGACACAGATATTTCTTCCACCTCCTTGGTGGGGCAGACATTTTTCTCTCTGT  
 GATGGCGTTTGACTGCTACATGGCCATCTCCAAGCCCTGCACTATGTGACCATCATGAGTAGAGGGCAA  
 TGCACCTGCCCTCATCTCTGCCTCTTGATGGGGGCTTTGTCCACTCCATCGTGAGATCTCCCTGTTGC  
 TGCCTCTCCCTTTCTGTGGACCAATGTTCTTGACACTTTCTACTGCGATGTCCCCAGGTCTCAAACCT  
 CACTTGCACTGACACTTTTGCTCTTGAGTTCTTGATGATTTCCAACAATGGCCTGGTCACTACCTGTGG  
 TTTATCTTCTGCTTGTGCTACACAGTCATCCTAATGACGCTGAGGTCTCAGGCAGGAGGGGGCAGGA  
 GGAAAGCCATCTCCACTTGACCTCCACATCACTGTGGTGACCTGCATTTTGTGCCCTGCATCTATGT  
 CTATGCCCGGCCCTTCACTGCCCTCCACAGAAAAGGCCATCTGTACCTTCACTGTCATCTCCCT  
 CTGCTGAACCTTTGATCTACACTCTGAGGAACAGGAAATGAAGTCAGCCATGAGAAGACTGAAGAGAA  
 GACTCGTGCCTTCTGAAAGGGAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Protein Sequence:</b> | <p>&gt;RC218114 representing NM_001004706</p> <p>Red=Cloning site Green=Tags(s)</p> <p>MELGNVTRVKEFIFLGLTQSQDQSLVFLFLCLVYMTLLGNLLIMVTVTCE SRLHTPMYFLLRNLAILD<br/> ICFSSTTAPKVLLDLLSKKKTISYTSMTQIFLHLLGGADIFSLVMAFDCYMAISKPLHYVTIMSRGQ<br/> CTALISASWMGGFVHSIVQISLLLPLFCGPNVLDTFYCDVPQVLKLTCTDTFALEFLMISNNGLVTTLW<br/> FIFLLVSYTVILMTLRSQAGGRRKAISTCTSHITVVTLHFVPCIVVYARPFTALPTEKAISVTFTVISP<br/> LLNPLIYTLRNQEMKSAMRRLLKRRLLVPSERE</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p> |
| <b>Species:</b>          | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Serotype:</b>         | AAV-2                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>ACCN:</b>             | NM_001004706                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>ORF Size:</b>         | 933 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Buffer:</b>           | PBS with 0.001% Pluronic F68                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Stability:</b>        | AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.   |
| <b>RefSeq:</b>           | <u>NM_001004706.1, NP_001004706.1</u>                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>RefSeq Size:</b>      | 936 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>RefSeq ORF:</b>       | 936 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Locus ID:</b>         | 219986                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>UniProt ID:</b>       | <u>Q8NGI4</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Cytogenetics:</b>     | 11q12.1                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>MW:</b>               | 34.8 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |