

## Product datasheet for **RC218163L1V**

### JMJD6 (NM\_001081461) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	JMJD6 (NM_001081461) Human Tagged ORF Clone Lentiviral Particle
Symbol:	JMJD6
Synonyms:	PSR; PTDSR; PTDSR1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001081461
ORF Size:	1242 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218163).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001081461.1</a>
RefSeq Size:	5445 bp
RefSeq ORF:	1245 bp
Locus ID:	23210
UniProt ID:	<a href="#">Q6NYC1</a>
Cytogenetics:	17q25.1
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS
MW:	47.4 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a nuclear protein with a JmjC domain. JmjC domain-containing proteins are predicted to function as protein hydroxylases or histone demethylases. This protein was first identified as a putative phosphatidylserine receptor involved in phagocytosis of apoptotic cells; however, subsequent studies have indicated that it does not directly function in the clearance of apoptotic cells, and questioned whether it is a true phosphatidylserine receptor. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]