

## Product datasheet for **RC218179L4V**

### STAU2 (NM\_014393) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | STAU2 (NM_014393) Human Tagged ORF Clone Lentiviral Particle   |
| Symbol:                   | STAU2  |
| Synonyms:                 | 39K2; 39K3   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-mGFP-P2A-Puro (PS100093)  |
| Tag:                      | mGFP   |
| ACCN:                     | NM_014393  |
| ORF Size:                 | 1437 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC218179).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_014393.1</a> , <a href="#">NP_055208.1</a>  |
| RefSeq Size:              | 4161 bp  |
| RefSeq ORF:               | 1440 bp  |
| Locus ID:                 | 27067  |
| UniProt ID:               | <a href="#">Q9NUL3</a>   |
| Cytogenetics:             | 8q21.11  |
| Domains:                  | DSRM   |
| Protein Families:         | Transcription Factors  |



[View online »](#)

**MW:** 52.8 kDa

**Gene Summary:** Staufen homolog 2 is a member of the family of double-stranded RNA (dsRNA)-binding proteins involved in the transport and/or localization of mRNAs to different subcellular compartments and/or organelles. These proteins are characterized by the presence of multiple dsRNA-binding domains which are required to bind RNAs having double-stranded secondary structures. Staufen homolog 2 shares 48.5% and 59.9% similarity with drosophila and human staufen, respectively. The exact function of Staufen homolog 2 is not known, but since it contains 3 copies of conserved dsRNA binding domain, it could be involved in double-stranded RNA binding events. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2009]