

Product datasheet for **RC218191L1V**

Heterogeneous Nuclear Ribonucleoprotein (A1 like) (HNRNPA1L2) (NM_001011724) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Heterogeneous Nuclear Ribonucleoprotein (A1 like) (HNRNPA1L2) (NM_001011724) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Heterogeneous Nuclear Ribonucleoprotein (A1 like)
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001011724
ORF Size:	960 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218191).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001011724.1 , NP_001011724.1
RefSeq Size:	2365 bp
RefSeq ORF:	963 bp
Locus ID:	144983
UniProt ID:	Q32P51
Cytogenetics:	13q14.3
Protein Pathways:	Spliceosome
MW:	34.2 kDa



[View online »](#)

Gene Summary:

Involved in the packaging of pre-mRNA into hnRNP particles, transport of poly(A) mRNA from the nucleus to the cytoplasm and may modulate splice site selection.[UniProtKB/Swiss-Prot Function]