

## Product datasheet for **RC218197L2V**

### Nestin (NES) (NM\_006617) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Nestin (NES) (NM_006617) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NES
Synonyms:	Nbla00170
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_006617
ORF Size:	4863 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218197).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006617.1</a>
RefSeq Size:	5591 bp
RefSeq ORF:	4866 bp
Locus ID:	10763
UniProt ID:	<a href="#">P48681</a>
Cytogenetics:	1q23.1
Protein Families:	ES Cell Differentiation/IPS
MW:	177.3 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the intermediate filament protein family and is expressed primarily in nerve cells. [provided by RefSeq, Sep 2011]