

## Product datasheet for RC218207L3V

### SLC12A8 (NM\_024628) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC12A8 (NM_024628) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC12A8
Synonyms:	CCC9
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_024628
ORF Size:	2142 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218207).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024628.4</a>
RefSeq Size:	3482 bp
RefSeq ORF:	2145 bp
Locus ID:	84561
UniProt ID:	<a href="#">A0AV02</a>
Cytogenetics:	3q21.2
Protein Families:	Druggable Genome, Transmembrane
MW:	78 kDa


[View online »](#)

**Gene Summary:**

This gene is thought to be a candidate for psoriasis susceptibility. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined. [provided by RefSeq, Sep 2010]