

Product datasheet for RC218219L2

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

HNF6 (ONECUT1) (NM_004498) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: HNF6 (ONECUT1) (NM_004498) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: HNF6

Synonyms: HNF-6; HNF6A

Mammalian Cell None

Selection:

Vector: pLenti-C-mGFP (PS100071)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC218219).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





 $[\]ensuremath{^*}$ The last codon before the Stop codon of the ORF.

ACCN: NM_004498

ORF Size: 1395 bp





HNF6 (ONECUT1) (NM_004498) Human Tagged Lenti ORF Clone - RC218219L2

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 004498.1</u>

RefSeq Size:1925 bpRefSeq ORF:1398 bpLocus ID:3175

UniProt ID: Q9UBC0
Cytogenetics: 15q21.3

Protein Families: ES Cell Differentiation/IPS, Transcription Factors

Protein Pathways: Maturity onset diabetes of the young

MW: 50.8 kDa


Gene Summary: This gene encodes a member of the Cut homeobox family of transcription factors. Expression

of the encoded protein is enriched in the liver, where it stimulates transcription of liver-expressed genes, and antagonizes glucocorticoid-stimulated gene transcription. This gene may influence a variety of cellular processes including glucose metabolism, cell cycle regulation, and it may also be associated with cancer. Alternative splicing results in multiple


transcript variants. [provided by RefSeq, Dec 2012]



Product images:



Circular map for RC218219L2



Double digestion of RC218219L2 using Sgfl and Mlul $\,$