

# **Product datasheet for RC218236**

### SOX1 (NM\_005986) Human Tagged ORF Clone

### **Product data:**

**Product Type:** Expression Plasmids

**Product Name:** SOX1 (NM\_005986) Human Tagged ORF Clone

Tag: Myc-DDK

**Symbol:** SOX1

Mammalian Cell Neomycin

Selection:

**Vector:** pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

ORF Nucleotide >RC218236 representing NM\_005986

Sequence: Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC

GCCGCGATCGCC

ATGTACAGCATGATGATGGAGACCGACCTGCACTCGCCCGGCGCCCCAGGCCCCCACGAACCTCTCGG CCAGGACCGGGTCAAACGGCCCATGAACGCCTTCATGGTGTGCCCCGCGGGCAGCGGCGCAAGATGGCC CAGGAGAACCCCAAGATGCACAACTCGGAGATCAGCAAGCGCCTGGGGGCCCGAGTGGAAGGTCATGTCCG TTACAAGTACCGGCCGCCGCAAGACCAAGACGCTGCTCAAGAAGGACAAGTACTCGCTGGCCGGCGGG CGGCCGTGGGCCAGCGCTGGAGAGCCCAGGCGGCGGCGGCGGCGGCTACGCGCACGTCAACGGCTG GGCCAACGGCGCCTACCCCGGCTCGGTGGCGGCGGCGGCGGCGCCGCGCCATGATGCAGGAGGCGCAG ACCACCGCACGCGCACCCGCACACCCGCAGCCCATGCACCGCTACGACATGGGCGCGCTGCAGTACAG CCCCATCTCCAACTCGCAGGGCTACATGAGCGCGTCGCCCTCGGGCTACGGCGGCCTCCCCTACGGCGCC CGGCCGCGTCGTCGGCCCCCTGGGCGCCTGGGCTCTCTGGTGAAGTCGGAGCCCAGCGGCAGCCCGCC CGCCCCAGCGCACTCGCGGGCGCCCTGCCCCGGGGACCTGCGCGAGATGATCAGCATGTACTTGCCCGCC ACTACCAGGGCGCGGGCGCGGGCGTGAACGGCACGGTGCCCCTGACGCACATC

AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA



**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Protein Sequence: >RC218236 representing NM\_005986

Red=Cloning site Green=Tags(s)

MYSMMMETDLHSPGGAQAPTNLSGPAGAGGGGGGGGGGGGGGGGGAKANQDRVKRPMNAFMVWSRGQRRKMA QENPKMHNSEISKRLGAEWKVMSEAEKRPFIDEAKRLRALHMKEHPDYKYRPRRKTKTLLKKDKYSLAGG LLAAGAGGGGAAVAMGVGVGVGAAAVGQRLESPGGAAGGGYAHVNGWANGAYPGSVAAAAAAAAAMMQEAQ LAYGQHPGAGGAHPHAHPAHPHHPHAHPHNPQPMHRYDMGALQYSPISNSQGYMSASPSGYGGLPYGA AAAAAAAAGGAHQNSAVAAAAAAAAASSGALGALGSLVKSEPSGSPPAPAHSRAPCPGDLREMISMYLPA GEGGDPAAAAAAAAQSRLHSLPQHYQGAGAGVNGTVPLTHI


**SGPTRTRRL**EQKLISEEDLAANDILDYKDDDDK**V** 

Chromatograms: <a href="https://cdn.origene.com/chromatograms/mg2396">https://cdn.origene.com/chromatograms/mg2396</a> e12.zip

**Restriction Sites:** Sgfl-Rsrll

**Cloning Scheme:** 





<sup>\*</sup> The last codon before the Stop codon of the ORF

**ACCN:** NM\_005986

ORF Size: 1173 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



#### **Reconstitution Method:**

- 1. Centrifuge at 5,000xg for 5min.
- 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
- 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** <u>NM 005986.3</u>

 RefSeq Size:
 4108 bp

 RefSeq ORF:
 1176 bp

 Locus ID:
 6656

 UniProt ID:
 000570

 Cytogenetics:
 13q34

Protein Families: Adult stem cells, ES Cell Differentiation/IPS, Induced pluripotent stem cells, Transmembrane

MW: 38.8 kDa

**Gene Summary:** This intronless gene encodes a member of the SOX (SRY-related HMG-box) family of


transcription factors involved in the regulation of embryonic development and in the

determination of the cell fate. The encoded protein may act as a transcriptional activator after

forming a protein complex with other proteins. In mice, a similar protein regulates the


gamma-crystallin genes and is essential for lens development. [provided by RefSeq, Jul 2008]


# **Product images:**



Circular map for RC218236







Western blot validation of overexpression lysate (Cat# [LY401812]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC218236 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).

Coomassie blue staining of purified SOX1 protein (Cat# [TP318236]). The protein was produced from HEK293T cells transfected with SOX1 cDNA clone (Cat# RC218236) using MegaTran 2.0 (Cat# [TT210002]).